

DATASHEET

An ELISA for the Detection of Antibodies to Trastuzumab in Human Serum

BACKGROUND

Trastuzumab (Herceptin®) is a humanized recombinant monoclonal antibody that selectively binds to the extracellular domain of the human epidermal growth factor receptor 2, HER2, a transmembrane protein overexpressed in 25-30% of breast cancers. Trastuzumab is used for the treatment of primary breast cancers which overexpress HER2. Herceptin binds with high affinity (kDa 5nM) to tumor cells over expressing HER2 resulting in loss of malignant growth and metastasis.

Evaluation of the comparability of trastuzumab biosimilars to the innovator drug should follow the guidelines laid out by the FDA and EMA. The analysis should be multifactorial, taking into account both the physicochemical characteristics and clinical performance of the biosimilar compared to the innovator. Eurofins Bioanalytical Services offers a full range of off-the-shelf trastuzumab assays for comparability testing of biosimilars including:

- PK assay
- ADA assay
- Nab assay
- Comparability testing
 - FcR & C1q binding
 - HER2 kinetic binding assay
 - ADCC assay


Trastuzumab has a mean half-life of 5.8 days in studies using a loading dose of 4 mg/kg dose followed by a weekly maintenance dose of 2 mg/kg. Mean trough and peak concentrations are approximately 79 µg/mL and 123 µg/mL respectively between weeks 16 and 32.

The purpose of this study was to develop a sensitive, specific and precise ELISA for the detection of antibodies to trastuzumab in human serum to support immunogenicity testing for clinical studies.

METHOD

The analytical method is an Affinity Capture Elution (ACE) ELISA where trastuzumab is used for capture of the anti-drug antibodies (ADA) in human sera. Unbound material is washed away and following acid dissociation and neutralization the captured ADA is transferred to a second plate and detected with biotinylated trastuzumab, followed by streptavidin-HRP, washed and visualized using TMB.

RESULTS



REFERENCES

1. Gopi Shankar, Viswanath Devanarayanb, Lakshmi Amaravadi, Yu Chen Barrett, Ronald Bowsher, Deborah Finco-Kentf, Michele Fiscella, Boris Gorovits, Susan Kirschner, Michael Moxnessj, Thomas Parishk, Valerie Quarmbyl, Holly Smith, Wendell Smith, Linda A. Zuckermano, Eugen Koren. Recommendations for the validation of immunoassays used for detection of host antibodies against biotechnology products. *Journal of Pharmaceutical and Biomedical Analysis* 48 (2008) 1267–1281.
2. James S. Bourdage, Carolyn A. Cook, Daphne L. Farrington, Jana S. Chain, Robert J. Konrad. An Affinity Capture Elution (ACE) assay for detection of anti-drug antibody to monoclonal antibody therapeutics in the presence of high levels of drug. *Journal of Immunological Methods* 327 (2007) 10–17.
3. Edith A. Perez, Edward H. Romond, Vera J. Suman, Jong-Hyeon Jeong, Nancy E. Davidson, Charles E. Geyer Jr, Silvana Martino, Eleftherios P. Mamounas, Peter A. Kaufman, and Norman Wolmark. *Journal of Clinical Oncology*. 29(25): 3366–73 (2011).
4. Tan, AR; Swain SM. "Ongoing adjuvant trials with trastuzumab in breast cancer". *Seminars in Oncology* 30 (5 Suppl 16) (2002): 54–64. doi:10.1053/j.seminoncol.2003.08.008.
5. Reflection Paper on Guidance for Laboratories that Perform the Analysis or Evaluation of Clinical Trial Samples. European Medicines Agency. EMA/INS/GCP/532137/2010
6. Guideline on bioanalytical method validation. EMEA/CHMP/EWP/192217/2009. Effective: February 2012.
7. Coussens L, Yang-Feng TL, Liao Y-C, Chen E, Gray A, McGrath J, et al. Tyrosine kinase receptor with extensive homology to EGF receptor shares chromosomal location with neu oncogene. *Science*. 1985; 230:1132–9.
8. Slamon DJ, Godolphin W, Jones LA, Holt JA, Wong SG, Keith DE, et al. Studies of the HER2/neu proto-oncogene in human breast and ovarian cancer. *Science*. 1989; 244:707–12.
9. Hudziak RM, Lewis GD, Winget M, Fendly BM, Shepard HM, Ullrich A. p185HER2 monoclonal antibody has antiproliferative effects in vitro and sensitizes human breast tumor cells to tumor necrosis factor. *Mol Cell Biol*. 1989; 9:1165–72.
10. Lewis GD, Figari I, Fendly B, Wong WL, Carter P, Gorman C, et al. Differential responses of human tumor cell lines to anti-p185HER2 monoclonal antibodies. *Cancer Immunol Immunother*. 1993; 37:255–63.
11. Baselga J, Norton L, Albanell J, Kim Y-M, Mendelsohn J. Recombinant humanized anti-HER2 antibody (Herceptin®) enhances the antitumor activity of paclitaxel and doxorubicin against HER2/neu over expressing human breast cancer xenografts. *Cancer Res*. 1998; 58:2825–31.
12. Similar biological medicinal products. CHMP/437/04 Sep 2005. CHMP/437/04 Rev. 1 May 2013.
13. European Medicines Agency. Similar biological medicinal products containing monoclonal antibodies: non-clinical and clinical issues. EMA/CHMP/BMWP/403543/2010. Jun 2012.
14. Similar biological medicinal products containing biotechnology-derived proteins as active substance: non-clinical and clinical issues. EMEA/CHMP/BMWP/42832/2005 Feb 2006. EMEA/CHMP/BMWP/42832/2005 Rev. 1 Jun 2013.
15. Similar biological medicinal products containing biotechnology-derived proteins as active substance: quality issues. EMEA/CHMP/BWP/49348/2005 Feb 2006. EMA/CHMP/BWP/247713/2012 Jun 2014.
16. Revision of the guideline on immunogenicity assessment of biotechnology-derived therapeutic proteins. EMA/275542/2013 Mar 2014.
17. United States FDA. Guidance for industry: Scientific Considerations in Demonstrating Biosimilarity to a Reference Product. Feb 2012.
18. United States FDA. Guidance for Industry: Quality Considerations in Demonstrating Biosimilarity to a Reference Protein Product. Feb 2012.
19. United States FDA. Guidance for Industry: Clinical Pharmacology Data to Support a Demonstration of Biosimilarity to a Reference Product. May 2014