

DATASHEET**An ELISA for the Detection of Antibodies to
Adalimumab in Human Serum****BACKGROUND**

Adalimumab (Humira[®]) is a humanized recombinant human IgG1 monoclonal antibody specific for tumor necrosis factor-alpha (TNF- α) that has been clinically applied for treatment of rheumatoid arthritis (RA), psoriatic arthritis, juvenile idiopathic arthritis, ankylosing spondylitis, chronic psoriasis, Crohn's disease, and ulcerative colitis.

Adalimumab mean steady-state trough concentrations of approximately 5 $\mu\text{g/mL}$ to 12 $\mu\text{g/mL}$ were observed at doses of 40 mg, 80 mg, and 160 mg of adalimumab. Among adalimumab monotherapy recipients, 3—26% had reported development of antibodies to adalimumab.

Evaluation of the comparability of adalimumab biosimilars to the innovator drug should follow the guidelines laid out by the FDA and EMA. The analysis should be multifactorial, taking into account both the physicochemical characteristics and clinical performance of the biosimilar compared to the innovator. Eurofins Bioanalytical Services offers an off-the-shelf ADA assay for comparability testing of biosimilar adalimumab.

The purpose of this study was to develop a sensitive and precise drug tolerant assay for the detection of antibodies to adalimumab in human serum. An Affinity Capture Elution (ACE) ELISA method was developed that is able to detect anti-adalimumab antibodies tolerant up to 20 $\mu\text{g/mL}$ adalimumab in human serum for supporting comparative analysis of biosimilar product clinical development, pharmacovigilance and therapeutic drug monitoring.

METHOD

The analytical method is an Affinity Capture Elution (ACE) ELISA where adalimumab is used for capture of the adalimumab anti-drug antibodies (ADA) present in human serum. Unbound material is washed away and following acid dissociation and neutralization, the captured ADA is transferred to a second plate and detected with biotinylated adalimumab, followed by streptavidin-HRP, washed and visualized using TMB.

RESULTS




FIGURE I. Representative curve for sensitivity of anti-adalimumab antibody in human serum.

TABLE I. Summary of assay performance results

Performance characteristic	Results
Sensitivity	4 ng/mL
Cut Point Assessment	30 individual normal male & female sera samples Floating screening cut point factor was established at 1.45
Selectivity (Matrix recovery)	9 out of 10 lots of normal human serum spiked with 250 ng/mL and 1250 ng/mL of positive control were within $\pm 25\%$ of reference
Precision Intra-assay Inter-assay	$\leq 13.6\% \text{ CV}$ $\leq 14.8\% \text{ CV}$
Drug Tolerance	250 ng/mL ADA is detectable in the presence of 20 $\mu\text{g/mL}$ of adalimumab

REFERENCES

1. Humira (adalimumab) package insert. North Chicago, IL: Abbott Laboratories; 2014 Sep.
2. Barrera P, Joosten LA, den Broeder AA, van de Putte LB, van Riel PL, van den Berg WB., Effects of treatment with a fully human anti-tumour necrosis factor alpha monoclonal antibody on the local and systemic homeostasis of interleukin 1 and TNFalpha in patients with rheumatoid arthritis. *Ann Rheum Dis.* 60 (2001) 660-669.
3. Taylor PC., Anti-tumor necrosis factor therapies. *Curr Opin Rheumatol.* 13 (2001) 164-169.
4. Karmiris K, Paintaud G, Noman M, Magdelaine-Beuzelin C, Ferrante M, Degenne D, Claes K, Coopman T, Van Schuerbeek N, Van Assche G, Vermeire S, Rutgeerts P. .Influence of trough serum levels and immunogenicity on long-term outcome of adalimumab therapy in Crohn's disease. *Gastroenterology.*137(2009):1628-1640.
5. Bartelds GM1, Wijbrandts CA, Nurmohamed MT, Stapel S, Lems WF, Aarden L, Dijkmans BA, Tak PP, Wolbink GJ. Clinical response to adalimumab: relationship to anti-adalimumab antibodies and serum adalimumab concentrations in rheumatoid arthritis. *Ann Rheum Dis.* 2007 Jul;66(7):921-6. Epub 2007 Feb 14.
6. Bartelds GM, Krieckaert CM, Nurmohamed MT et al. Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. *JAMA* 2011; 305(14):1460-68.
7. Gopi Shankar, Viswanath Devanarayanb , Lakshmi Amaravadi , Yu Chen Barrett, Ronald Bowsher, Deborah Finco-Kentf, Michele Fiscella, Boris Gorovits, Susan Kirschner, Michael Moxnessj, Thomas Parishk, Valerie Quarmbyl, Holly Smithm,Wendell Smithn, Linda A. Zuckermano, Eugen Koren. Recommendations for the validation of immunoassays used for detection of host antibodies against biotechnology products. *Journal of Pharmaceutical and Biomedical Analysis* 48 (2008) 1267–1281.
8. James S. Bourdage, Carolyn A. Cook, Daphne L. Farrington, Jana S. Chain, Robert J. Konrad. An Affinity Capture Elution (ACE) assay for detection of anti-drug antibody to monoclonal antibody therapeutics in the presence of high levels of drug. *Journal of Immunological Methods* 327 (2007) 10–17.
9. Reflection Paper on Guidance for Laboratories that Perform the Analysis or Evaluation of Clinical Trial Samples. European Medicines Agency. EMA/INS/GCP/532137/2010
10. Guideline on bioanalytical method validation. EMEA/CHMP/EWP/192217/2009. Effective: February 2012.